

γ H2AX Foci Counting Workflow

Author: Albert Wang

Last updated: 2025-07-14

Summary


This document describes a semi-automated MATLAB pipeline for quantifying γ H2AX foci within individual nuclei. The workflow covers image import, preprocessing, segmentation, foci detection, per-cell metric computation, and results export.

Workflow Steps

MCF10A_CLR_1h_1 is used as example image here


1. Image Import

- Use `uigetfile` to select two TIFF files: first the nucleus channel, then the γ H2AX (foci) channel. Read both images with `imread`
- Convert each RGB image into a single-channel grayscale intensity map by summing the red, green, and blue channels (`I_nucleus`, `I_foci`).



2. Thresholding and Mask Creation

- Prompt the user for:
 - A nucleus threshold value (0–255, used with `imbinarize`).
 - A foci threshold (0–255 intensity units).
 - Thresholding mode:
 - Dynamic (1): Determines a per-nucleus foci threshold.
 - Universal (2): Uses a single global foci threshold.
- Nucleus mask:
 - Binarize `I_nucleus` using `imbinarize`.
 - Clean the binary mask with `bwareaopen` and `imfill`.
 - Use `bwconncomp` to identify connected nuclei.
 - Use `regionprops` to extract nucleus area and centroid.
 - Discard nuclei smaller than a defined area threshold (e.g., 1,000 pixels) and record remaining indices for downstream analysis.
- Foci mask:
 - Create an initial foci mask by thresholding `I_foci`.
 - For each retained nucleus:
 - Apply the selected thresholding mode to define foci regions.
 - Remove small specks (< 2 px) with `bwareaopen`.
 - Aggregate per-nucleus foci masks into a global foci mask.
 - For each discarded (small) nucleus:
 - Remove surrounding foci in a square region (default: 90 px half-width) to eliminate false positives.
- Can also use MATLAB's Color Thresholder App for auto thresholding (<https://www.mathworks.com/help/images/image-segmentation-using-the-color-threshold-app.html>)



3. Per-Nucleus Metric Computation


For each retained nucleus:

- Extract the grayscale intensity values from `I_foci`.

- Compute summary statistics:
 - Mean, median, mode, standard deviation.
- Apply per-nucleus foci thresholding and count foci:
 - Use bwconncomp and regionprops to exclude small objects.
- Store results as a row of:
 - Index, Area, MeanInt, MedianInt, ModeInt, StdInt, FociCount

4. Results Export

- Store all per-nucleus results in a MATLAB table.
- Write the table to an Excel file using writetable.
- Generate a 2x2 segmentation figure including:
 - Nucleus segmentation overlay
 - Foci segmentation overlay
 - Color-coded nucleus index map (with overlaid numbers)
 - Original foci image
- Save the figure as a high-resolution JPEG (print -djpeg -r600).



Note: Cells with their nuclei extremely close to each other (CHLA-20 in particular) would be difficult for the program to distinguish individual nucleus (see example below). Discard these data to avoid skewing the overall trend.

